

ProKlenz™ Booster and ProKlenz NpH Sterile Detergents

Concentrated detergents for controlled environments

The ProKlenz Booster and ProKlenz NpH Sterile Detergents are designed to remove organic and inorganic residues from product and non-product contact surfaces. Sterile filtered and filled in an ISO Class 5 aseptic environment, each package is double bagged for easy entry into your controlled environments and aseptic core areas. These products have an extensive documentation package to meet your validation objectives and are supported by the STERIS Technical Services team for individualized application and validation assistance.




Features	Benefits
ProKlenz Booster Sterile Detergent contains surfactants and hydrogen peroxide	Utilizes multiple cleaning mechanisms to clean a broad spectrum of soils
ProKlenz NpH Sterile Detergent contains multiple surfactants	Effectively cleans soils with a mild, formulated surfactant system
Moderate pH	Can be used in manual cleaning and compatible with many substrates
Double bagged	Allows easy entry into controlled environments
0.2 micron filtered, filled in ISO Class 5 (Grade A, Grade B at rest) environment and each lot is USP <71> sterility tested	Introduces minimal bioburden into your facility
Manufactured under a FDA cGMP, ISO 9001 and ISO 13485 compliant quality system	Supports quality assurance and quality control
REACH and European Union Regulation (EC) No. 648/2004 on Surfactants compliant ingredients	An excellent choice for global harmonization of cleaning processes.

Applications

- Residue removal from:
 - Product contact surfaces
 - Floors
 - Walls
 - Ceilings
 - Countertops
 - BSL hoods
 - Carts and other equipment
- Products may be applied by:
 - Mopping
 - Wiping



Life Sciences



An important aspect to any microbial control program is a strategy to address disinfectant and sporicidal residues. Since any disinfectant or sporicide will leave some residue, we recommend that every program include a rinsing or residue removal step. For facilities that do not routinely rinse or have tough residues, we recommend using a detergent. Our ProKlenz Booster and ProKlenz NpH Sterile Detergents are ideal for this application because they can be easily brought into controlled environments and are formulated to remove challenging soils. Our Technical Service group can help customize your contamination control program to take advantage of this excellent product.

Ordering Information

ProKlenz Booster Sterile Detergent

World wide

2.6 ounce	24 per case	1S1337WR
1 gallon	4 per case	1S1308WR

ProKlenz NpH Sterile Detergent

World wide

2.6 ounce	24 per case	1S1237WR
1 gallon	4 per case	1S1208WR

About STERIS Corporation

STERIS Corporation is a leading provider of infection prevention and surgical products and services, focused primarily on critical healthcare, pharmaceutical and research markets around the world. The Company supplies a broad array of equipment, consumable and service solutions that help assure productivity and quality. The Company is listed on the New York Stock Exchange under the symbol STE. For more information, visit www.steris.com.

STERIS has a comprehensive offering of detergents, disinfectants, skin care products and sterility assurance products that support your needs. We also have world class technical support to design the most effective cleaning program for your facility.

STERIS Corporation
5960 Heisley Road
Mentor, OH 44060
USA

T +1 800 444 9009

STERIS Offices Worldwide

Asia Pacific	+65 (0) 68 41 7677
Canada	+ 1 800 444 9009
France	+33 (0)2 38 84 85 40
Germany	+49 (0) 221 466120
Italy	+39 (0) 2 21303 424
Latin America	+1 800 444 9009
Spain	+34 (0) 916 585 920
United Kingdom	+44 (0) 1256 840400

www.STERISLifeSciences.com

© 2014 STERIS Corporation
All rights reserved

No part of this publication may be reproduced in any material form without the written permission of STERIS Corporation.

LSFC-SS4088.A-EN-E May 2014